

## Introducing the new West-Call i-Dome®



### Overview

West-Call's® new *Room Status* hardware is a standalone device that displays patient and room status — not just bed status. *Room Status* is about the size of a light switch panel, and is installed just outside of a patient's door. It has one clear LED light that can display a variety of colors (for example, purple as shown above), as well as three configurable buttons below it. *Room Status* is used to alert medical personnel of patient status via its LED for easy viewing from a hallway. *Room Status* has a totally customizable array of light assignments spanning a vast assortment of colors. It also operates independently of any nurse call system, providing flexibility and customization well beyond existing, proprietary solutions. *Room Status* is fed from external systems and does not actively monitor, protecting the privacy of the patient. Light color assignments are non-standard; protecting the privacy rights of patients for whom *Room Status* provides alerts. *Room Status* is controlled primarily through a highly customizable web interface using an IP-based system that runs on CAT5. **Multiple simultaneous states that cycle colors can also be used providing flexibility and robust optimization for any healthcare environment.**

### Key Features

- Can be configured to display patient or room status.
- Can show status from other devices such as beds.
- Entirely stand alone (not dependent upon a nurse call system).
- Data is received from external systems, does not actively monitor.
- Fully customizable web interface.
- Uses an IP-based system running on CAT5 cable.
- Create any color to light up the LED light.
- Supports multiple states; can cycle between the colors to show multiple status issues.