



EV NMS

The Endeavor Nurse Call System comes complete with a 15" CRT, an IBM compatible computer, keyboard, mouse, and a phone. An optional 15" touch screen flat panel display is also available (shown above). The interactive interface uses a call list, displaying active calls and staff location. Large descriptive buttons with colorful animated graphics make finding and using functions simple. Help buttons located throughout the Endeavor program allow the user to receive answers to operational questions on demand. Full duplex voice communication is available from the master station to patient stations, staff stations, duty stations, and other master stations. Other features include fully integrated radio paging, staffing reports, staff assignment, and patient information, all customizable to meet the nursing unit's needs.

- ▶ Full Duplex Speech
- ▶ Integrated Radio Paging (optional)
- ▶ Customizable Features
- ▶ Online Help Menus
- ▶ Interactive Call List
- ▶ Staff Reports
- ▶ Staff Assignments
- ▶ Patient Information
- ▶ Color 15" CRT or optional 15" Touch Screen Flat Panel Display

DESCRIPTION	MODEL NO.
Endeavor Master Station with 15" Touch Screen LCD, Phone, Keyboard, Mouse, and Computer	EV NMS/LCD
Endeavor Master Station with 15" CRT, Phone, Keyboard, Mouse, and Computer	EV NMS

OPTIONS	
TQI for call reporting Feature	/TQI
Battery Back Up and Surge Protection	UPS-8500

- To combine features, specify corresponding code letters that begin with a (/) to the end of the model number (e.g., EV NMS/TQI)

MONITOR

Standard Monitor	15" CRT
Optional Monitor	15" LCD Touch Screen (Five-Wire Resistive)
Power Consumption	30W
Resolution	1024 x 768 at 60, 70 or 74 Hz
Response Time	30ms
Brightness	300cd/m ²
Dimensions (with base)	Width: 15.4", Height (Inc. Base): 15.4", Depth (Inc Base): 8.2"
Weight	Actual: 15.2 lb (6.9 kg) Shipping: 21.6 lb (9.8 kg)
Touch Screen	AccuTouch Five-Wire Resistive

COMPUTER

Power Supply	External DC Power Adapter
Motherboard	All-In-One Board
Processor	1.8GHz Intel processor, 128K Cache
Memory	128 MB SDRAM
Video Card	Integrated. Intel Extreme Graphics
Connectivity	One 9-Pin D-Sub Serial Port, One 25-Pin D-Sub Parallel Port, One RJ-45 Network Connector, PS/2 Mouse & Keyboard Port, 15-Pin VGA Port, 4 USB on back 2 USB on front
Hard Drive	20GB, EIDE, 9.5 MM, 5400 RPM
CD-ROM Drive	Removable CD-ROM Drive
Floppy Disk Drive	Optional
Voltage	90-135 V at 60 Hz; 180-265 V at 50 Hz
Weight	7.2 lbs.
OTHER	
Phone	DC 12V
Keyboard	84 Keys
Cabling	Category 5 Cable From Master Computer To Endeavor Central. West-Call® Special Black Cable From Master Phone to Central



View West-Com's complete West-Call® line of healthcare communications solutions at www.westcall.com

West-Com Nurse Call Systems Incorporated
P.O. Box 640 • Danville, CA 94526 • 800-761-1180 • fax: 925-831-0842 • info@westcall.com