

Now Available

Crest-brand pillow speakers from West-Com Nurse Call Systems



West-Call Pillow Speakers are built on Crest Healthcare Supply's award-winning Excel platform and are now available to all West-Call distributors. Call West-Com to purchase your demo sample today and present a superior quality product to your customers.

West-Call Pillow Speakers are available in Digital, Analog and Infrared styles. All styles feature a range of standard functions, allowing you to choose the configuration that meets the needs of your customers.



- Drop tested onto concrete for extreme durability
- Bottom entry cable with strain relief for fewer repairs
- Award-winning features and design
- Superior sound quality that lasts
- Ergonomic design for easy handling

Call West-Com today to place your order!



Infrared

The only pillow speaker that combines complete TV, DVD and VCR control.

Digital (10 Key & Big Button)

Designed for use with digital TVs.

Analog

Designed for use with non-digital TVs.

OTHER QUALITY PRODUCTS ALSO AVAILABLE

- CALL CORDS
- SPECIALITY NURSE CALL
- ADAPTERS
- RECEPTACLE KITS

- WALL PLATES
- JUMPER CABLES
- BED COMMUNICATION CABLES



Call today for great values on quality nurse call products

phone: 800.761.1180 or 707.428.5900
fax: 707.428.5941

West-Com Nurse Call Systems, Inc.
2200 Cordelia Road, Fairfield, CA 94534



2200 Cordelia Road Tel (800) 761-1180
Fairfield, CA 94534 Fax (707) 428-5941

Introducing A Wireless Sound Solution!

The new IR Audio System breaks down the barriers of wiring to give you the ability to use pillow speakers in new places!



Patent Pending

IR Audio™ System




New IR Audio™ System

Experience high-quality pillow speaker sound in any healthcare environment without the wiring! The new IR Audio System gives you the ability to use West-Com's proprietary EntertainIR™ Pillow Speaker without the need for wiring between the TV and the nurse call station.

Pillow speakers have been used for many years in hospital patient rooms to provide nurse call, lighting management, TV control and a private TV listening

experience. They help reduce excessive environmental noise from TV use by placing the speaker within close proximity of the user. They can also condense the number of cords and controls each person needs into a single, convenient device. However, many long term care facilities and some areas in hospitals have been unable to use pillow speakers because they lack the wiring required in the walls and ceiling to carry the sound from the TV to the nurse call station.

- **Make pillow speaker use available in any room without wiring**
- **No pillow speaker interface needed on the TV**
- **Same great features and sound quality as Excel pillow speakers**
- **Sets up in minutes**
- **Works with many nurse call systems**



The new IR Audio System eliminates this barrier by carrying sound signals wirelessly from the TV to a Receiver which connects to your existing nurse call station. The proprietary IR Audio System eliminates the line of sight requirements typical of IR communication systems to ensure that uninterrupted crisp, clear sound is delivered to the EntertainIR™ Pillow Speaker at all times regardless of room traffic and won't interfere with other signals.

Quick, Easy Set-Up

Setting up the IR Audio System is easy using three basic components. The Transmitter mounts near the TV and is plugged into your TV using the headphone, line out or pillow speaker audio output and a power source. The Receiver is mounted on a wall near the nurse call station and connects directly to the station and a power source. The EntertainIR™ Pillow Speaker is then plugged into the Receiver to deliver sound and to control the TV functions with IR signals just like your home remote.

Call West-Com at **1-800-761-1180** to learn more about the IR Audio™ System.



2200 Cordelia Road Tel (800) 761-1180
 Fairfield, CA 94534 Fax (707) 428-5941

Award Winning Excel EntertainIR™ Pillow Speaker

now controls your cable or satellite TV set top boxes!



- > Video on demand control option
- > Menu navigation control
- > DVD control option

Innovative design provides improved capabilities and flexibility that other manufacturers can't offer.

West-Call® Excel EntertainIR™ Pillow Speakers **Now with Set Top Box control**

West-Com now offers pillow speakers with the capability to operate cable or satellite TV Set Top Boxes. This allows you to give your patients access to these devices as well as video on demand features and menu navigation, without the need for a separate remote control. Set Top Box pillow speakers are built on our award-winning EntertainIR™ platform and carry all the great features and benefits of our Excel family. If you're currently using a pillow speaker and a remote control to operate your TV and Set Top Box, call West-Com to find out more about how you can get a pillow speaker that does it all!

- **Drop tested onto concrete for extreme durability**
- **Superior sound quality that lasts**
- **Bottom-entry cable with strain relief for fewer repairs**
- **Award-winning features and design**



Call West-Com at **1-800-761-1180** to learn more about set top box control with Excel EntertainIR™.