

**EV-SPS1**

The West-Call® Endeavor Patient Station is a single-gang device that allows patients to place a routine call and control the television and room lights through a pillow speaker. It also works with a standard call cord. When the cord is removed from the station, it will place a call. The red button on the station lights when a call is placed and cancels the call when pressed.

Features

- ▶ Places a routine call
- ▶ Works with pillow speakers or call cords
- ▶ Places a call when cord is removed from station
- ▶ Large call cancel button
- ▶ Fits into a one gang back box
- ▶ Cord-out override function
- ▶ Capable of controlling the room's TV and room lights when used with the appropriate pillow speaker
- ▶ Call assurance light on station
- ▶ Places a call if station loses power

Ordering Information

DESCRIPTION	MODEL NO.
West-Call® Endeavor Single Patient Station with 8-Pin modular jack	EV-SPS1
PILLOW SPEAKERS & CALL CORDS	
Pillow Speaker with 8' Cable, TV, NC, L1, L2	31010-003
Pillow Speaker with 8' Cable, TV, NC	31010-004
Pillow Speaker with 8' Cable, TV, NC, L1, L2, Headset Jack	31010-101
Nurse Call Cord with 8' Cable, 10-P Mod Jack	CC96-661
Single Breath Call Cord with 10 Pin Connector	GEB-175

- For use with the West-Call® Endeavor system

Technical Specification

Operating Voltage	12-30 VDC; 12-18 VAC
Power Consumption	1 watt maximum
Dimensions	4 ½"H x 2 ¾" W x 1 ½" D
Mounting	4-inch square backbox with one-gang plaster ring, Raco 190 with 768 one-gang adapter or equivalent. Minimum 2 1/8" depth required. Normal mounting is vertical. For remodel areas, use Raco #500 or 503 with #977 Grip-Lok's or equivalent. Minimum 2 1/2" depth required
Faceplate	Anodized aluminum
Lamp Life	10,000 hour minimum
Switch Life	100,000 cycles minimum
Cable Requirement	Use 6 conductor #22AWG cable
Termination Requirement	Plug-on connector allows rapid mass termination. Mating connector is a Panduit CE100F22-5
Certifications	UL-1069 for US and Canada, E139082

NOTE: This information should not be used for final design. Consult West-Com for total system requirements. Specifications subject to change without notice.



View West-Com's complete West-Call® line of healthcare communications solutions at www.westcall.com

West-Com Nurse Call Systems, Inc.
2200 Cordelia Road • Fairfield, CA 94534 • 800-761-1180 • info@westcall.com