



Features:

- Programmable call annunciation; calls can be filtered from any Controller and/or Central.
- Provides visual call display and sound only; audio handset not included.
- Software runs on durable IPC-3 computer.
- 15-inch touch screen monitor allows for easy call/list management (Patients & Staff).
- Configurable call priority determines the order of call display.
- Active calls can be muted, while any new incoming calls will produce sound.
- Mute button is 'activated' only when a call is received.
- Capable of interfacing with paging system(s) to send alert pages to staff members.
- Uses the same .wav files for call tones and call priority ordering as a FocusCare® master station.
- Patient information can be provided from Data Manager via ADT feed or can be manually entered.
- Staff information can be sent from another master station or entered manually.
- Use with the West-Call® Novus®, Odyssey with FocusCare® Ver. 1.1.8 & Endeavor systems.
- Connects via Ethernet to the Nurse Call network.
- Compatible *only* with FocusCare® Version 1.1.8 platform.



Computer (IPC-3) Specifications*

Input Voltage	12-24VDC; 12VDC @ 3A, 24VDC @ 1.5A
Power Consumption	Typical: 24W, Maximum: 36W
Motherboard	Integrated Intel 855GME
Processor	Low Voltage 1.4GHz Intel Pentium®-M
Memory	512 MB RAM
Storage	2GB Solid-state Flash
Dimensions	10.41" x 2.72" x 5.4" (264.5x69.2x137.25mm)
Weight	5 LBS
Operating Temperature	-20°~60°C (-4°~140°F), RH: 95% @ 40°C non-condensing

Monitor Specifications*

Display Size	15" Diagonal
Touch Screen	IntelliTouch Surface Wave technology with durable, scratch-resistance glass screen
Resolution	1280 x 1024 max
Response Time	18 ms (typical)
Brightness	280-nit brightness
Power Management	Meets TCO '03 and Energy Star® standards
Power Supply	10.41" x 2.72" x 5.4" (264.5x69.2x137.25mm)
AC input (line) voltage	90-265VAC, 50/60Hz, 50W maximum
Dimensions	Width: 15.2", Height: 14.4", Depth: 2.6"
Environment (Operating)	0°~40°C (32°-104°F), RH: 20% - 80% non-condensing
Warranty	3 years

Accessories

Mouse	USB Optical
Keyboard	Mini USB
Cabling	Ethernet and AC Power x 2

**Note: This information should not be used for final design. Consult West-Com for total system requirements. Specifications subject to change without notice.*



View West-Com's complete West-Call® line of healthcare communications solutions at www.westcall.com

West-Com Nurse Call Systems, Inc.
2200 Cordelia Road • Fairfield, CA 94534 • 800-761-1180 • info@westcall.com