



LIAISON™

Bridging All Your Communication Needs...

LIAISON™ FEATURES

Liaison™ is a general purpose flexible interface between West-Com's nurse call products and external products, services, and data.

Liaison™ not only provides interfaces to paging systems, two-way wireless phones (such as SpectraLink, Ascum, Vocera and Voalte) but can also provide connections between different data sources and targets. For example, a data stream from a TCP/IP connection to a telemetry monitoring system can be selected as a source, with the output being a Novus® Nurse Call System to create "calls" based on input trigger parameters specified for the input stream.



REVOLUTIONIZING
the patient care environment

WC WEST-CALL

West-Com Nurse Call Systems, Inc.

2200 Cordelia Road, Fairfield, CA 94534 1.800.761.1180
info@westcall.com www.westcall.com

Inputs:

West-Call® Nurse Call IP protocol • TCP/IP • TAP • Binary or Text formatted records • Record definition tool (i.e., assign field names to fields in variable or fixed records)

Output Protocols Supported:

(Individual licensing for each protocol):

TAP • ESPA • SpectraLink OAI • SNPP • SMTP
Vocera VMI • Voalte • Cisco Nurse Connect
Fixed Record • Delimited Record • Aeroscout • SQL

Note: This product replaces the need for the following products: Radio Page Server (RPS) and Telephone Line Server (TLS).

- ✓ Synchronizes with FocusCare® DataManager for interface/device definitions.
- ✓ Multiple transports selectable (IP or Serial)

* Some features are planned but not currently implemented. Always work with West-Com to properly design and install your system. Features and specifications are subject to change without notice.

