



**Nurse Aide  
Locator**



**Yellow and Red  
Level Locator**



**Green and Yellow  
Level Locator**

The Locator station provides two levels of staff location, in addition to other unique features such as: cancellation of calls and selective audio page functions. When used with the West-Call® Odyssey Plus or West-Call® Endeavor systems, locators can be a signaling device for any purpose.

- ▶ Highly visible call LED to assure the call has been placed; 500,000 hour life expectancy
- ▶ Easily operated push color-coded area
- ▶ Fire retardant high impact molded thermoplastic
- ▶ Solid state circuitry
- ▶ Plug-in connectors for easy system connection/service
- ▶ Designed to work with West-Call® Odyssey and Endeavor
- ▶ Custom colors and text available

### VERSIONS AND MODEL NUMBERS:

STATION TYPE	ODYSSEY	ENDEAVOR
Two Level Locator with a Green Level 1 Button and a Yellow Level 2 Button	LOC-1000	EV LOC
Two button locator with a Green Nurse Button and a Yellow Aide Button	LOC-1000/NA	EV LOC/NA
Two Level Locator with a Yellow Level 1 Button and a Red Level 2 Button	LOC-1001	EV LOC/1

- For use with the West-Call® Odyssey Plus and West-Call® Endeavor systems
- Custom colors and text available upon request

Operating Voltage	28 VDC
Power Consumption	1 watt maximum
Size	4 1/2"H x 2 3/4"W x 1 1/4"D
Mounting	4" square backbox with 1 gang plaster ring, Raco 190 with 768 or equivalent. 2 1/8" minimum depth required. Normal mounting is <u>Vertical</u>
Faceplate	Hospital grade ABS UL94V-0 molded thermoplastic with inlaid membrane
Faceplate Color	Warm Grey (Pantone 1C)
Membrane	.010 Lexan with hardcoat finish
Membrane Color	.010 Lexan with hardcoat finish background color Grey (Pantone 4C). Copy color Warm Grey (Pantone 1C and 11C)
Termination	Plug on connector allows rapid mass termination. Mating connector - Panduit CE156F22-10
Cable Requirements	6 conductor #22AWG cable to associated station or I/O board
Certifications	UL-1069 for US and Canada, E139082

**NOTE:** This information should not be used for final design. Consult West-Com for total system requirements. Specifications subject to change without notice.



View West-Com's complete West-Call® line of healthcare communications solutions at [www.westcall.com](http://www.westcall.com)

West-Com Nurse Call Systems, Inc.  
2200 Cordelia Road • Fairfield, CA 94534 • 800-761-1180 • [info@westcall.com](mailto:info@westcall.com)