



### NCM-FC/H

The West-Call® Nurse Call Master Station comes complete with an IPC-3 Computer, 15" flat panel touch screen display with integrated phone, mini keyboard, mouse, Windows® XP Embedded operating system and FocusCare® License.

- Interactive touch screen user interface with a graphical floor map displaying all rooms assigned to the nursing unit.
- Large descriptive buttons with colorful animated graphics allow for simple identification and utilization of functions.
- Full duplex voice communication from the master station to patient stations, staff stations, staff duty stations and other master stations.
- Full messaging, staffing reports, staff assignment and patient information.
- Highly customizable to meet the nursing unit's needs.

- ▶ Full Duplex Audio
- ▶ 15" Flat Panel Touch Screen Display
- ▶ Interactive Room Map Display
- ▶ Staff Reports
- ▶ Master to Master Communication
- ▶ Staff Assignments
- ▶ Patient Information
- ▶ Customizable Features
- ▶ Integrated Messaging

**VERSIONS AND MODEL NUMBERS:**

STATION TYPE	MODEL NO.
West-Call® Nurse Call Master. IPC-3, 15" LCD touch panel with integrated phone, keyboard, mouse and FocusCare® License.	NCM-FC/H
OPTIONS	
Extended Warranty (5-year total)	IPC-3W

- For use with the West-Call® Novus®, Odyssey and Endeavor Systems
- Contact a West-Com factory representative for additional information on West-Call® Master Stations

**MONITOR**

Display Size	15" Diagonal
Power Management	Meets TCO'03 and ENERGY STAR® standards
Resolution	1280 x 1024 max
Response Time	18 ms (typical)
Brightness	280-nit brightness
Dimensions	Width: 15.2", Height: 14.4", Depth: 2.6"
Touch Screen	IntelliTouch Surface Wave technology with durable, scratch-resistant glass screen
Power Supply	AC input (line) voltage: 90-265 VAC, 50/60 Hz; 50 W max.
Environment	0°C to 40°C; RH:
Warranty	3 years

**ACCESSORIES**

Computer	IPC-3
Input	Mini keyboard, optical mouse
Phone	Integrated handset or stand-alone
Software	West-Call® FocusCare® and support software
Warranty	2 years

**OTHER**

Cabling	Ethernet, AC power x2, and West-Call® Special Black
Certifications	UL-1069 for US and Canada, E139082

**NOTE:** This information should not be used for final design. Consult West-Com for total system requirements. Specifications subject to change without notice.



View West-Com's complete West-Call® line of healthcare communications solutions at [www.westcall.com](http://www.westcall.com)

West-Com Nurse Call Systems, Inc.  
2200 Cordelia Road • Fairfield, CA 94534 • 800-761-1180 • [info@westcall.com](mailto:info@westcall.com)