



P-Tag

The West-Call® LPS P-Tag has been designed to be worn by patients, using feedback from hospital personnel. The unique construction consists of a reusable electronic core and a disposable, single use (and waterproof) outer shell which makes the P-Tag safe, comfortable and economical to use.

The P-Tag clam shell snaps easily on and safely attaches to any standard patient wristband, or similar band or lanyard.

- ▶ Attached to regular patient wristband
- ▶ Batteries: 120,000 transmissions
- ▶ Reusable Core
- ▶ Single patient use outer shell
- ▶ Disposable outer shell with batteries
- ▶ Waterproof (shower, bath, etc.)

VERSIONS AND MODEL NUMBERS:

DESCRIPTION	MODEL NUMBER
P-Tag Core. Reusable electronics for Patient Tag	30175003
P-Tag Shell Kit. Disposable shell with batteries for Patient Tag	30175004

- For use with the West-Call® Novus®, Odyssey and Endeavor systems

Operating Voltage	3.0 VDC, 140 mAh @ 20 mA
Dimensions	Core 1.15 x 0.54 in Shell 1.37 x 1.45 x 0.69 in.
Material	K-Resin®, KR01
Material Color	Light grey (Pantone 427C)
Weight	Total 0.437 oz.
Certifications	UL-1069 for US and Canada, E139082

NOTE: This information should not be used for final design. Consult West-Com for total system requirements.
Specifications subject to change without notice.



View West-Com's complete West-Call® line of healthcare communications solutions at www.westcall.com

West-Com Nurse Call Systems, Inc.
2200 Cordelia Road • Fairfield, CA 94534 • 800-761-1180 • info@westcall.com