

NOVUS[®] VIP



Rated Number ONE Nurse Calling System

West-Call's Novus[®] received the highest user experience overall ranking. Tested criteria included reliability, service and support, functionality, interfacing and ease of use.

We Make Nurse Call Easy

Easiest to install.
Easiest to learn.
Easiest to operate.
Easiest to service.

Most Reliable System Made in the USA

West-Call products are made in the USA. Always have been. Always will be. Backed with a five-year warranty. West-Call products' reliability surpasses our competition.



West-Com Nurse Call Systems, Inc.

NOVUS® i-BED Status Annunciator

West-Com's patented i-DOME® is installed in the VIP annunciator panel to give you a "bed status" panel.



i-DOME®

i-Status®

- i-Dome® lens indicates a change in bed status.
- Room lens indicates type of bed status change and the room in which the change has occurred.
- i-Dome® displays i-Bed status from the West-Call® Novus® Visual bed connector.
- CAT 5 Cabling.



Pull Cord

Visual Patient Station

Pull Cord, Push for Help



Dome Light

West-Com Introduces

NOVUS® VIP

Our New Smart Visual System



- The highly visible i-Dome® lens offers "unlimited" color combinations to indicate the level of call (White=Routine Call, Flashing White=Bathroom Emergency Call, Red=Emergency Call, Flashing Blue=Code Blue).
- 32 Room identification lenses that are customizable to each facility's needs.

- Room level calls also light according to the level of call.
- The system comes with 16 standard tones. However, custom file tones are available upon request.
- Speaker ■ Power light ■ CAT 5 Cabling



Code Blue

Code Pink

Staff Emergency

Staff Assist

Code White

Staff Emergency/
Code Blue

NOVUS® Remote Mini-Master

The Novus® VIP panel can be installed in the hallway to inform staff of patient room status or to direct staff to the appropriate patient room when used as a Mini-Master, or in a staff area when used as an expanded duty station in conjunction with a full Novus® system. The VIP panel can also be installed as a "stand alone" annunciator in a visual only system.



Novus® Patient Station with True Duplex Audio

The Mini-Master indicates:

- Type of call displayed on i-DOME® lens, (White=Routine Call, Flashing White=Room Emergency, Red=Staff Emergency, Flashing Yellow=Bed Exit, Flashing Blue=Code Blue) and many more.
- Room lens shows location of call and color of lens light indicates type of call.
- Power indication LED.
- Speaker grill (comes with 16 standard tones).
- CAT 5 Cabling.
- This mini-master works in conjunction with the full Novus® master station, indicating types of call associated with each room.



Novus Master Station