



West-Com's Novus[®] + Gen4 Pillow Speakers =
**EXTRAORDINARY
COMBINATION**

NOVUS® + GEN 4 = *An Extraordinary Pair!*

We think our new Novus® Nurse Call System is pretty extraordinary. It features a new computer, Cat-5 cabling, the new i-Dome® for instant status of at-risk patients, and newly-designed patient and staff stations.

Our new Gen4 pillow speakers are just as extraordinary. They were designed to operate all the advanced features of the latest Nurse Call systems, like the Novus®, while being easy to operate and comfortable to hold. Gen4 combines the best features of previous pillow speakers, suggestions from customers, and new materials and techniques, making it a step ahead of anything else in the nurse call industry today.

- **Infection Control**

Our proprietary design helps prevent infectious contaminants from getting into the pillow speaker, where they could be passed from patient to patient. In addition, Gen4's smooth


shape allows for quick and easy cleaning, ensuring that dirt and grime cannot collect on the unit.

- **Patient Satisfaction**

Gen4's case is streamlined and lightweight, and we made sure the speaker provided crystal clear sound. All the buttons are easy to find and easy to press, making it as comforting to use as it is comfortable to hold.

- **Oxygen Safety**

Don't be fooled by pillow speakers that are labeled "UL Recognized." The UL Recognized mark does not ensure the product has undergone oxygen safety testing or that the nurse call accessory has been tested with the nurse call system it is intended to be used with. Gen4 Pillow Speakers are UL Listed and tested with West-Call Nurse Call Systems and verified for safe use in an oxygen-enriched atmosphere. For more information on oxygen safety, visit www.oxygensafe.com



NOVUS® + GEN 4 = *Flexible, Efficient and Powerful Solutions*

Direct or Digital

Gen4 Pillow Speakers are available in Digital and Direct Channel Access models for all West-Call Nurse Call Systems.

No "Cord Out"

West-Com's exclusive breakaway cable gives added protection to the patient station and eliminates false "cord out" calls.

Optional Stereo Headphone Jack

Ideal for creating a quieter environment for patient rooms.

Reprogrammability

Gen4 pillow speakers can be easily reprogrammed in the field for different televisions.

Optional Analog Volume Control

Ideal for situations where two pillow speakers need to control one television.

Contact your West-Call dealer to find out how to replace your current nurse call equipment with something truly extraordinary!



GEN 4 Features

Large Buttons

Gen4's large buttons are recognizable at a glance and easy to press. The simplified layout makes contacting the Nurse and controlling TV functions simple.

Direct Request™ Buttons

Pain Med and Bed Pan buttons let you take advantage of the nurse call system's advanced features.

Volume Control Options

All models come standard with digital volume controls, and are available with analog volume for situations where two pillow speakers need to control one TV.

Flex And Strain Relief

Gen4's long flex relief was designed for longer life and greater

flexibility. The smooth design facilitates easy cleaning and resists dirt and grime buildup. Internal strain relief prevents damage to internal components.

Easy-to-Clean Case

The streamlined design facilitates quick and easy cleaning.

Reprogrammable for Different TVs

There's no need to manage numerous pillow speakers to support various hospital TV brands in your facility. Each Gen4 Pillow Speaker can be quickly reprogrammed.

Auxiliary Functions

The Gen4 platform allows for up to six auxiliary buttons to control lighting, temperature, drapes, nurse call features, etc.



REVOLUTIONIZING *the patient care environment*



West-Com Nurse Call Systems, Inc. is dedicated to providing the healthcare industry with communications solutions that improve the quality of patient care. The Gen4 Pillow Speaker, paired with Novus® Nurse Call System, enable facilities to improve productivity while enhancing the patient experience. For more details and a demo, please contact us at 800.761.1180. Visit West-Com's full line of healthcare communications solutions for medical offices, long-term care facilities and acute-care facilities at www.westcall.com



West-Com Nurse Call Systems, Inc.

2200 Cordelia Road, Fairfield, CA 94534 1.800.761.1180
info@westcall.com www.westcall.com