

Bed	Admit Dt/Tm	P.I.	Staff	Gender	Age
2208-A	09/22	JE	ND SA	M	37

Bed	Admit Dt/Tm	P.I.	Gender	Age	Diagnosis	NPO	Acuity	Status	
2208-A	09/22	JE	ND SA	M	37	Thrombosis	N	4	Routine
2208-B	09/24	HR	ND SA	M	5 M	Coma	N	3	Routine
2211-A	09/25	VF	MM SA	F			N	5	Emergency
2211-B	09/21	ML	MM SA	F	46	Strep	N	5	Routine
2212-A	09/24	BS	MM SA	F			N	5	Routine
2214-A	09/22	RE	KW CC BA	M	27	Heart Attack	Y	5	Routine
2215-A	09/22	AG	KW CC BA	M	26	Broken Leg	N	2	Routine
2218-A	09/22	MC	SB KH	M	28		N	5	Routine
2218-B	09/24	JH	SB KH	M	33		N	1	Routine
2219-A	09/25	CG	SB KH	F	1 W		N	5	Routine
2220-A	09/22	AS	BB JD KH	F	46		N	2	Routine
2220-B	09/24	BS	BB JD KH	F		Flu	N	5	Routine
2221-A	09/21 05:03	PB	BB JD KH	M	8		N	5	Routine
2221-B	09/25	RS	BB JD KH	M	86		N	5	Routine

iBed State

Bed	Admit Dt/Tm	P.I.	State	Gender	Age	NPO	Status	Acuity	Staff
2208-A	09/22	JE	[Exit]	M	37	N	Routine	4	ND SA
2208-B	09/24	HR	[Exit]	M	5 M	N	Routine	3	ND SA
2211-A	09/25	VF	[Exit]	F		N	Emergency	5	MM SA
2211-B	09/21	ML	[Exit]	F	46	N	Routine	5	MM SA
2212-A	09/24	BS	[Exit]	F		N	Routine	5	MM SA
2214-A	09/22	RE	[Exit]	M	27	N	Routine	5	CC KW BA
2215-A	09/22	AG	[Exit]	M	26	N	Routine	2	CC KW BA
2218-A	09/22	MC	[Exit]	M	28	N	Routine	5	SB KH
2218-B	09/24	JH	[Exit]	M	33	N	Routine	1	SB KH
2219-A	09/25	CG	[Exit]	F	1 W	N	Routine	5	SB KH
2220-A	09/22	AS	[Exit]	F	46	N	Priority	5	BB JD KH
2220-B	09/24	BS	[Exit]	F		N	Routine	5	BB JD KH
2221-A	09/21 05:03	PB	[Exit]	M	8	N	Priority	3	BB JD KH
2221-B	09/25	RS	[Exit]	M	86	N	Routine	5	BB JD KH

White Board

The FocusCare® White Board displays current status information for a set of beds. It is intended for use in hallways or nursing areas where a large screen monitor can provide current assignments, patient status information and various other custom selectable fields. This functionality can replace the traditional 'grease boards' used in many facilities. Any changes to the bed assignments, staffing, patient info, bed availability, etc. made through the web interface, FocusCare® masters, ADT or other methods will automatically be reflected on the White Board display.

With the optional Stryker® iBed™ interface, the FocusCare® White Board is capable of displaying data directly from each bed such as bed exit, brake, side rail, fowler, and bed height alarms.

- ▶ Viewable on high resolution wall-mounted displays (not included)
- ▶ Customize which data is actually displayed
- ▶ Customize the colors, fonts and views
- ▶ Bed status and patient information are easily located on screens throughout the patient floor areas

VERSIONS AND MODEL NUMBERS:

DESCRIPTION	MODEL NO.
<i>White Board</i> interface installed on an IPC-3 computer. Monitor/Display sold separately	White Board

- For use with the West-Call® Novus®, Odyssey and Endeavor systems (Odyssey only when using FocusCare® software)
- For minimum PC requirements, please refer to IPC-3 data sheet