



 Check EOL Directive  Isolation  Allergy  Bed Set  Fall Risk

Lighting The Way To Superior Care

i-DOME[®]

i-DOME *Takes Room Status to a Higher Level*

Similar in size to a light switch, once this small device is installed outside a patient's room, it can be programmed to alert staff to a variety of urgent conditions through its interface with the nurse call system, patient bed, and wireless devices. Its clear LED light and three configurable buttons can display a vast assortment of colors in a customizable array of light assignments. West-Com designed the i-DOME® to offer a multiplicity of options for improving patient care - from pain management and documentation to reduction of falls - by alerting staff to patients identified as "at-risk". Nurses can

activate timers to visually remind them when it's time to reassess a patient's medication. Staff are instantly alerted to a change in bed status. Housekeeping efficiency and productivity increases too, as i-DOME® notifies staff that it's time to clean a patient's room, and Admissions gets notification once the room is ready. An excellent tool to support compliance with Joint Commission guidelines (JACHO), i-DOME® offers facilities a simple, yet powerful way to improve patient care in the medical environment.



i-Bed Awareness & **i-DOME**

- Allow you to focus more on patient care by reducing time spent on routine tasks such as bed checks.
- Visually remind you to reset the Chaperone® Bed Exit System when the patient returns to bed.
- Notify staff of any changes in the *i*Bed Awareness settings.

Room Status is fed from external systems and does not actively monitor, protecting the privacy of the patient. Light color assignments are non-standard, protecting the privacy rights of patients for whom Room Status provides alerts. Room Status is controlled primarily through a highly customizable web interface using an IP-based system that runs on CAT5.



i-DOME *Features*

- Can be configured to display patient or room status.
- Indicates status from other devices such as beds.
- Entirely stand-alone (not dependent upon a nurse call system).
- Data is received from external systems; does not actively monitor.
- Fully customizable web interface.

- Is an IP-based system running on CAT5 cable.
- Select colors red, green, blue, purple and yellow to light up the LED indicator light.
- Supports multiple states; can cycle between the colors to show multiple status levels.



Stryker S3™ Bed



REVOLUTIONIZING[®] *the patient care environment*



West-Com Nurse Call Systems, Inc. is dedicated to providing the healthcare industry with communications solutions that improve the quality of patient care. West-Com's i-DOME takes room status hardware to a higher, more sophisticated level, offering facilities a simple, yet powerful way to improve patient care and the medical environment. For more details and a demo, please contact us at 800.761.1180. Visit West-Com's full line of healthcare communications solutions for medical offices, long-term care facilities and acute-care facilities at www.westcall.com



West-Com Nurse Call Systems, Inc.

2200 Cordelia Road, Fairfield, CA 94534 1.800.761.1180
info@westcall.com www.westcall.com

