



### i-Status®

West-Call's® new *i-Status*® is a stand-alone device that displays patient and room status in a manner that is clear to staff while maintaining patient privacy. It can be used for such things as Pain Medication Assessment reminders, Clean/Dirty Bed status, and Assistance calls. The *i-Status*® device can virtually take the place of any of the current calls or tasks in the nurse call system. It operates independently of any nurse call system, patient station, or dome light, providing flexibility and customization well beyond existing, proprietary solutions. Light color assignments are user-defined, protecting the privacy rights of patients for whom *i-Status*® provides alerts. *i-Status*® is controlled primarily through a highly customizable web interface using an IP-based system that runs on CAT5, and it integrates seamlessly with the FocusCare® software. Multiple simultaneous states provide flexibility and robust optimization for any healthcare environment.

- ▶ Provides Status for Pain Medications, Bed Status (Clean/Dirty/In Progress), can virtually take the place of any call within a nurse call system
- ▶ CAT5 cabling
- ▶ Up to 63 devices within one network
- ▶ Can be independent of any Nurse Call System
- ▶ IP-based
- ▶ Customizable Web Interface

### VERSIONS AND MODEL NUMBERS:

STATION TYPE	NOVUS®
<i>i-Status</i> ® is controlled primarily through a highly customizable web interface using an IP-based system that runs on CAT5, and it integrates seamlessly with the FocusCare® software	<i>i-Status</i> ®

- Can interface with the West-Call® Novus® system only

Operating Voltage	28 VDC
Power Consumption	.06 Watts idle, .75 Watts maximum
Dimensions	4 1/2" H x 2 3/4" W x 2" D
Mounting	4" Square back box with 1 gang plaster ring, Raco 190 with 768 or equivalent. 2-1/8" minimum depth required. Normal Mounting is vertical.
Faceplate Material and Color	Hospital Grade Polylac PA765NAT with Light Grey Opaque
Membrane	0.008" Textured Polyester with MDWC-Brown 400
Cable Requirements	One CAT5 daisy chained through the i-Dome®/i-Status® network up to 63 devices.
Termination	Two RJ-45 connectors
Certifications	UL-1069 for US and Canada, E139082

**NOTE:** This information should not be used for final design. Consult West-Com for total system requirements. Specifications subject to change without notice.



View West-Com's complete West-Call® line of healthcare communications solutions at [www.westcall.com](http://www.westcall.com)

West-Com Nurse Call Systems, Inc.  
2200 Cordelia Road • Fairfield, CA 94534 • 800-761-1180 • [info@westcall.com](mailto:info@westcall.com)